


Applications of DNA sequencing

Overview

- DNA Sequencing Applications
- Reference Genome availability
- Quality Assessment and Trimming of DNA sequence reads

Whole Genome Shotgun Sequencing

- Start with genomic DNA
- DNA is sheared into fragments
 - Physical
 - Acoustic shearing (Covaris)
 - Sonication (Bioruptor)
 - Hydrodynamic force (Hydroshear)
 - Enzymatic (transposase, DNase I)
 - Chemical
- Ideally, would like a very uniform size selection
 - paired end: depends on kit, from 200-600bp
 - mate pairs: 3-20 Kbp



Many Uses


1. Resequencing analysis

We know a reference genome, and want to find *variants* (blue) in a background of errors (red)




2. Counting

We have a reference genome (or gene set) and want to know how *much* we have. Think gene expression/microarrays.



3. Assembly

We don't have a genome or any reference, and we want to construct one.
(This is how all new genomes are sequenced.)




The basic approaches to using sequencing data.

So Many Uses

Biologists are clever and have hijacked shotgun sequencing for many types of new biological assays, mostly having to do with clever ways of treating the DNA and preparing the libraries.

- RNASeq
 - convert RNA to cDNA for sequencing
 - we'll cover this more later
- Targeted DNASeq
 - only sequence regions of interest
- ChIP Sequencing
 - Sequence genome-wide binding sites
- Reduced complexity sequencing (GBS/RADSeq) for genotyping
 - sequence a smaller part of the genome
- Methylation sequencing/ bisulfite seq
 - Identify DNA bases that have methylation
- Hi-C
 - 3D Chromatin structure
- Etc.



Coverage (Depth of Coverage)

- All applications depend on appropriate coverage
- Depth of Coverage = the average number of reads that align to or “cover” a reference base
- Mean mapped read depth =

$$\frac{\text{\# of mapped bases}}{\text{\# of reference bases}}$$

- Often displayed as a histogram
- Coverage uniformity matters too

Which SNP is believable?

- Just by eye (there are real statistics to this that we'll learn more about later)

Ref : GAAGTGGCATATGGCTGTGAAGAAAAAG

R1 : GAACTGGCATATGGCTGTGA

R2 : ATATGGCAGTGAAGAAA

Which SNP is believable?

- Just by eye (there are real statistics to this that we'll learn more about later)

Ref : GAAGTGGCATATGGCTGTGAAGAAAAAG

R1 : GAACTGGCATATGGCAGTGA


R2 : AACTGGCATATGGCAG

R3 : ATATGGCAGTGAAGAAA

R4 : ATGGCAGTGAAGAAAA

R5 : GCAGTGAAGAAAAAG

Depth of Coverage as a Histogram



What do you think causes the long tail?

<https://www.biostars.org/p/48633/>

App 1: Whole Genome Resequencing

- Reference genome is already available
- Sequencing one or more individuals, usually from the same species
- Discover variations in the genomes between and within samples
 - Straightforward
 - SNPs
 - insertions
 - deletions
 - Harder
 - rearrangements
 - copy number changes

1. Resequencing analysis

We know a reference genome, and want to find *variants* (blue) in a background of errors (red)



Quality Assessment



Trimming



Quality Assessment



Mapping to a Reference



Visualization



Calling variants



Assessing functional impact of variants



Submit to SRA

Example
Workflow
(Pipeline)

App 2: GBS/RADSeq

Polymorphism discovery

Polymorphisms are an essential genomic tool for:

- Population Structure
 - Association mapping
 - Pedigree mapping
 - QTL mapping
 - Phylogeny
- Use the high volume and low cost of sequencing to replace SNP chips and microsatellites
 - Whole genome resequencing is expensive!
 - How to efficiently use NGS for discovering markers?
 - How to efficiently use NGS to do the genotyping?

Restriction Site-associated DNA sequencing (RADSeq)


- Random genome reduction
- Developed Baird et al 2008
- Identify and score thousands of genetic markers
- Subsampling only at specific sites defined by restriction enzyme
- Randomly distributed across the target genome
- From many individuals using Illumina technology


Patented by University of Oregon, licensed by:



Non profit use at universities can be licensed for free

Baird NA, Etter PD, Atwood TS, Currey MC, Shiver AL, et al. (2008) Rapid SNP Discovery and Genetic Mapping Using Sequenced RAD Markers. PLoS ONE 3(10): e3376.





1° reads	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTTG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTTG			
2° reads	TGC	G	GGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTTG	
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTC	C	TG	
	TGCAGGACACACAGGAGC	A	GAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTTG	
	TGCAGGACACACAGGAGCTGAGCCATTTC	T	TGCGGGCTCC	A	GACCAAACGTTTG	
	TGCAGGACA	G	ACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCA	G

ddRAD


- May have the problem of too many sites across the genome (even if you use a rare cutter)
- Need a way to more accurately control the number of loci sequenced
- Double digest RAD
- Peterson et al 2012
- Also patented

ddRAD

- Double digest RAD
- (Peterson et al 2012)
- Simpler and cheaper library construction
- restriction digest with two enzymes simultaneously
- eliminate random shearing and end repair
- explicitly use size selection
- Sequence fragments generated by cuts with both REs and which fall within the size-selection window


A

RAD sequencing




B

double digest RADseq



Peterson et al (2012) Double Digest RADseq: An Inexpensive Method for De Novo SNP Discovery and Genotyping in Model and Non-Model Species. PLoS ONE 7(5): e37135.

Its important to carefully select the technology that meets your goals.




Peterson et al (2012) Double Digest RADseq: An Inexpensive Method for De Novo SNP Discovery and Genotyping in Model and Non-Model Species. PLoS ONE 7(5): e37135.

GBS: Genotyping by Sequencing

- Elshire et al 2011
- Increased efficiency and cost benefits
- Reduced sample handling
- Methylation-sensitive REs used to filter out the repetitive portion of the genome
- Better barcoding system
- Fewer steps
- Free to use and to sell (ie not licensed)

*In 2015 I paid \$46 per sample
96-plex, 1 plate
Including SNP calling
Including enzyme optimization
External rate at Cornell



Workflow


- Informatics can be difficult if you don't have a genome
- Much easier if you have a reference genome – can be similar to whole genome resequencing
- If not, some software packages:
 - Universal Network Enabled Analysis Kit (UNEAK) – part of the TASSEL software suite
 - Stacks

 Stacks



App 3: Bisulfite Sequencing

- DNA methylation
 - First discovered epigenetic mark
 - methyl groups are added to DNA
 - Suppresses transcription when present in promoter
 - Adenine and Cytosine can be methylated in prokaryotes
 - Only cytosine is methylated in eukaryotes



Bisulfite treatment

- How to figure out where methylation occurs while sequencing?
- Treatment of DNA with bisulphite:
 - Unmethylated cytosine -> uracil
 - 5-methylcytosine stays the same
- Sequencing can yield single- nucleotide resolution of methylation patterns

Problems:

- Incomplete conversion
- DNA degradation during conversion
- 5-methylcytosine and 5-hydroxymethylcytosine both read as a C in bisulphite sequencing



Bisulfite conversion



PCR



Workflow

- Need special software for mapping
 - Bismark
 - BSMap
 - BSMapper
- Downstream analysis
 - MethyKit – statistics, visualization, tiling windows



Overview

- Many types of omics yield different types of information
- Your data analysis protocols must be tailored to fit the data production protocols

Reference Genomes

All the methods we talked about so far depend on (or are easier with) a reference genome.

How many genomes are out there to use as a base for mapping reads?

UCSC Genome Browser

- Mammal (55)
- Other Vertebrate (27)
- Deuterostome (3)
- Insect (13)
- Nematode (6)
- Other (2) – yeast and sea hare
- Viruses (1) - Ebola


UCSC Genome Browser on Ebola virus Sierra Leone 2014 (G3683/KM034562.1/eboVir3) Assembly


move <<< << < > >> >>> zoom in 1.5x 3x 10x base zoom out 1.5x 3x 10x 100x

KM034562v1:1-18,957 18,957 bp.

[hg38 replaces hg19 as default human assembly](#)

KM034562v1 KM034562v1 NP VP35 VP40 KM034562v1 GP VP30 VP24 KM034562v1 L





Reference Genomes

NCBI Genome

NCBI Genome has 4 levels:

- Complete - all chromosomes are gapless and have no runs of 10 or more ambiguous bases (Ns), there are no unplaced or unlocalized scaffolds, and all the expected chromosomes are present
- Chromosome - there is sequence for one or more chromosomes, gaps OK.
- Scaffold - some sequence contigs have been connected across gaps to create scaffolds, but the scaffolds are all unplaced or unlocalized
- Contig - nothing is assembled beyond the level of sequence contigs

NCBI Genome Records

- Viruses
 - 19,743 genomes
 - 19,084 complete (> 96%)
- Prokaryotes
 - 162,676 genomes
 - 11,576 complete (7.0%)
 - 2,057 chromosome level (1.3%)
- Eukaryotes
 - 6,563 genomes
 - 47 complete (< 1%)
 - 874 chromosome level (13.3%)

What are we doing today?

1. Examining Read Quality
2. Quality/Adapter Trimming

Quality Control

Goals

- Is my data of sufficient quality to use?
- The instrument assigns a confidence value to each base. Are the bases high quality overall?
- Does the complexity look normal?
 - PCR and library prep problems can lead to duplication of the same sequences over and over
- Are there adapters or other over-represented sequences?
- Are there lane batch effects?

FASTQC



Babraham Bioinformatics

Accepts input formats:

- FastQ (all quality encoding variants)
- GZip compressed FastQ
- SAM
- BAM

Does a 12 point analysis of quality



Generates an html output file

- ✓ [Basic Statistics](#)
- ✓ [Per base sequence quality](#)
- ✓ [Per tile sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ! [Per base sequence content](#)
- ✓ [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ✓ [Overrepresented sequences](#)
- ✓ [Adapter Content](#)
- ! [Kmer Content](#)

FASTQC



Babraham Bioinformatics




<http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>

Trimming

- From the quality control step, we know where the problems are
- All illumina reads tend to have degrading quality at the end of the read
- Get rid of the bad data, keep the good data
 - Cut adapter sequences from the read.
 - Trim off low quality bases
 - Drop a read entirely if is too low quality or too short

Trimmomatic

- Optimized for Illumina NGS
- Very flexible
- Handles paired end data well
- Threaded
- Detects adapter read through
- No read through:
- Read through:



Bolger, A. M., Lohse, M., & Usadel, B. (2014). Trimmomatic: A flexible trimmer for Illumina Sequence Data. *Bioinformatics*, btu170.


Trimmomatic

The current trimming steps are:

- **ILLUMINACLIP**: Cut adapter and other illumina-specific sequences from the read. Must specify adapter sequence. Comes with basic Illumina adapter files, make sure yours are in there or add yours!
- **SLIDINGWINDOW**: Perform a sliding window trimming, cutting once the average quality within the window falls below a threshold.
- **LEADING**: Cut bases off the start of a read, if below a threshold quality
- **TRAILING**: Cut bases off the end of a read, if below a threshold quality
- **CROP**: Cut the read to a specified length
- **HEADCROP**: Cut the specified number of bases from the start of the read
- **MINLEN**: Drop the read if it is below a specified length after trimming

Trimmomatic

- Maximum Information Quality Filtering:
- Retain low-quality bases early in a read in order to make sure the read is sufficiently long to be informative
- Trimming process becomes more strict later in the read



Skewer

- Faster than trimmomatic
- A bit less flexible
- “Gentle” quality trimming by default

<https://github.com/relipmoc/skewer>

Jiang, H., Lei, R., Ding, S.W. and Zhu, S. (2014) Skewer: a fast and accurate adapter trimmer for next-generation sequencing paired-end reads. BMC Bioinformatics, 15, 182.

Trimming

- Current community wisdom:
 - Quality trimming reduces error
 - But also reduces content and contiguity
- Gentle trimming is preferred – many times the defaults are too stringent, you will lose lots of data!
- Application matters
 - For mapping and counting, gentle to no trimming (phred 3 to 5)
 - For assembly and variant calling, a bit more trimming is good (phred 10 to 15)